


DDL Inc provides mechanical testing services and consulting to medical device manufacturers and the packaging industry. They own multiple TestResources systems shown in the applications below:


Contact information as follows
www.testedandproven.com

DDL Inc
10200 Valley View Road, Suite 101
Eden Prairie, Minnesota 55344
tel: 952.941.9226
fax: 952.941.9318
1.800.229.4235




**Catheter Body
Elongation Tensile
Test**


Tested with
TestResources 650M
Test System equipped
with G240K Vise Action
Grips




Syringe Test tested using a
TestResources 200R Test System



TestResources **G94**
Lightweight Air Grips
performing tensile tests
on thin film.



Surgical Blade Test
using a TestResources
100Q Test System




TestResources **200R Test System** mounted horizontally to perform needle puncture tests on bottles used in medical applications.



TestResources **G7-1 Roller Capstan Grips** used for lap shear tests on plastics.



TestResources **200R Test System with Collet Grips** perform Catheter Luer Fitting Tensile Strength Test




TestResources **G7W Wedge Action Grips** perform low force fatigue tests on welds and solder joints.



TestResources **G238 Bend Fixtures** testing ferrite rods.



TestResources **G240K Vise action grips** performing strength tests on stainless steel foil.



TestResources **G76-1 Fine Wire Capstan Grips** perform strength tests on small diameter wires common to medical devices.